



# Instruction Manual

## Hydraulic Hand Dynamometer with Digital LCD Gauge




Manufacturer and Master  
Distributor of Physical Therapy  
and Rehabilitation Products

**REF** 12-0247

ver 0723

## Parts/Specifications



<b>Specifications</b>	300 lb / 136 kg capacity
<b>Grip adjust range</b>	1.35 in / 3.35 cm
<b>Weight</b>	22.6 oz / 638 gm



## **Hydraulic Hand Dynamometer with Digital LCD Gauge**

### **Purpose / Intended Use**

The purpose of this Medical Device is to measure a patient's grip strength. It can be used to evaluate initial and ongoing impacts of hand trauma and dysfunction. The Baseline Hydraulic Hand Dynamometer can be used to measure isometric force and peak strength with five adjustable grip positions. This helps track strength and rehabilitation over time and is ideal for physical therapy professionals.

### **How to Use**

Set handle to comfortable grip for patient. Reset max indicator to zero. Have patient squeeze with maximum force, note reading. Reset to zero for next test.

### **Components**

- Machined aluminum handle, post and body
- Bronze bellows
- Stainless steel hydraulic tubing
- Teflon bushings
- Non-toxic mechanical hydraulic pump fluid
- Digital Gauge


Accuracy greater than 97%.

### **Data**


The Baseline® Hydraulic Hand Dynamometer with Digital LCD Gauge can utilize some of data pertaining to the Jamar® Hydraulic Hand Dynamometer. The internal workings of both are hydraulic and bellows-operated.

# Instructions

## How to turn on/off unit:




## How to zero out unit:




## How to view max:



## How to clear max:



## How to switch between LB and KG:



## Norms for Adult Grip Strength

A recent study by Dr. Virgil Mathiowetz indicates that "... individuals using the Baseline® dynamometer are justified in using the normative data collected with the Jamar® dynamometer..."

For each test of grip strength, the subject was seated with shoulder adducted and neutrally rotated with the elbow between 0° and 15° ulnar deviation.

The standard test protocol used the mean of three strength tests as a resultant score. A score was taken with both the dominant (right) and non-dominant (left) hands.

The rest results show a relationship between:


- hand strength vs. age
- hand strength of men vs. hand strength of women
- dominant hand strength vs. non-dominant hand strength

<b>Average Performance of all Subjects on Grip Strength (pounds)- Test results</b>											
Men							Women				
Mean	SD	SE	Low	High	Age	Hand	Mean	SD	SE	Low	High
121.0	20.6	3.8	91	167	20-24	R	70.4	14.5	2.8	46	95
104.5	21.8	4.0	71	150		L	61.0	13.1	2.6	33	88
120.8	23.0	4.4	78	158	25-29	R	74.5	13.9	2.7	48	97
110.5	16.2	4.4	77	139		L	63.5	12.2	2.4	48	97
121.8	22.4	4.3	70	170	30-34	R	78.7	19.2	3.8	46	137
110.4	21.7	4.2	64	145		L	68.0	17.7	3.5	36	115
119.7	24.0	4.8	76	176	35-39	R	74.1	10.8	2.2	50	99
112.9	21.7	4.2	73	157		L	66.3	11.7	2.3	49	91
116.8	20.7	4.1	84	165	40-44	R	70.4	13.5	2.4	38	103
112.8	18.7	3.7	73	157		L	62.3	13.8	2.5	35	94
109.9	23.0	4.3	65	155	45-49	R	62.2	15.1	3.0	39	100
100.8	22.8	4.3	58	160		L	56.0	12.7	2.1	37	83
113.6	18.1	3.6	79	151	50-54	R	65.8	11.6	2.3	38	87
101.9	17.0	3.4	70	143		L	57.3	10.7	2.1	35	76
101.1	26.7	5.8	59	154	55-59	R	57.3	12.5	2.5	33	86
83.2	23.4	5.1	43	128		L	47.3	11.9	2.4	31	76
89.7	20.4	4.2	51	137	60-64	R	55.1	10.1	2.0	37	77
76.8	20.3	4.1	27	116		L	45.7	10.1	2.0	29	66
91.1	20.6	4.0	56	131	65-69	R	49.6	9.7	1.8	35	74
76.8	19.8	3.8	43	117		L	41.0	8.2	1.5	29	63
75.3	21.5	4.2	32	108	70-74	R	49.6	11.7	2.2	33	78
64.8	18.1	3.7	32	93		L	41.5	10.2	1.9	23	67
65.7	21.0	4.2	40	135	75+	R	42.6	11.0	2.2	25	65
55.0	17.0	3.4	31	119		L	37.6	8.9	1.7	24	61
104.3	28.3	1.6	32	176	All Subjects	R	62.8	17.0	0.96	25	137
93.1	27.6	1.6	27	160		L	53.9	15.7	0.88	23	115


### References:

1. Gill D, Reddon J, Renney C, Stefanyk W: Hand Dynamometer: Effects of Trials and Sessions. *Perceptual and Motor Skills* 61: 195-8, 1985.
2. Everett P, Sills F: The relationship of Grip Strength to Stature, Somatotype Components, and Anthropometric Measurements of the Hand. *The Research Quarterly* 23: 161-6, 1952
3. Mathiowetz V, Federman S, Wiermer D: Grip and Pinch Strength: Norms for 6 to 19 Year Olds. *The American Journal of Occupational Therapy* 40: 705-11, 1986
4. Mathiowetz V, Donahoe L, Renells C: Effect of Elbow Position on Grip and Key Pinch Strength. *The Journal of Hand Surgery* 10A: 694-7, 1985
5. Mathiowetz V, Kashman N, Volland G, Weber K, Dove M, Rogers S: Grip and Pinch Strength: Normative Data for Adults. *Archives of Physical Medicine and Rehabilitation* 66: 69-74, 1985.

### Grip Strength (Dominant Hand)



### Grip Strength Difference Between Dominant and Non-Dominant Hand



\* charts generated from data published in Mathiowetz's article "Grip and Pinch Strength: Normative Data for Adults", Archives of Physical Medicine and Rehabilitation 66: 69-74, 1985

THE COLLEGE OF  
**ST. CATHERINE**  
2004 Randolph Avenue  
St. Paul, Minnesota 55105  
(612) 690-6000 FAX(612) 690-6024

Virgil Mathiowetz, PhD, OTR  
Associate Professor  
Department of Occupational Therapy  
College of St. Catherine  
2004 Randolph Avenue  
St. Paul, MN 55105-1794

March 18, 1993

Mr. Elliott Goldberg,  
Marketing Director  
Fabrication Enterprises Inc.  
Trent Building  
South Buckout Street  
Irvington, NY 10533

Dear Mr. Goldberg,

Recently, I completed the study to determine whether the Baseline and Jamar hydraulic dynamometers can be used interchangeable. A draft of the report has been completed and sent to you. In the summary, I concluded that, "The data from this study suggest that the Jamar and Baseline hydraulic hand dynamometers measure equivalently for practical purposes. As a result, individuals using the Baseline Dynamometer are justified in using the normative data, which was collected with the Jamar dynamometer (Mathiowetz et al., 1985; 1986)." This conclusion assumes that the same standard procedures are followed as were used in the original normative data studies.

Sincerely,



Virgil Mathiowetz, PhD, OTR  
Associate Professor &  
Research Consultant

## **WARRANTY**

The Baseline® Hydraulic Hand Dynamometer with Digital LCD Gauge is warranted for 2 years covering parts and labor from date of purchase. If unit needs repair, contact your local dealer or Fabrication Enterprises, Inc.

250 Clearbrook rd, Elmsford NY, 10523 (USA)

tel: 914-345-9300 • 800-431-2830

fax: 914-345-9800 • 800-634-5370



Fabrication Enterprises Inc.  
250 Clearbrook Rd, Suite 240  
Elmsford, NY 10523 (USA)  
tel: +1-914-345-9300 • 800-431-2830  
fax: +1-914-345-9800 • 800-634-5370  
FabEnt.com



AJW Technology Consulting GmbH  
Breite Strasse 3  
40213 Düsseldorf (Germany)